



### Defibrillator Maintenance Checklist

<b>Defibtech DDU-100</b>
Defibtech DDU-100 Serial Number: .....
Defibtech DDU-100 Location:.....

Date: Monthly check required. Refer to Defibtech User Manual kept with defibrillator.	Month: .....	Month: .....	Month: .....	Month: .....	Month: .....	Month: .....
Active Status Indicator(ASI) is flashing green						
Unit and accessories are free from for damage, dirt and contamination.  Clean or replace if necessary (refer to manual).						
Battery pack and pads have not past expiration dates						
Comments:						
<b>Checked by: (signature + print name )</b>						

**Copy of completed form to be sent to Occupational Nurse Consultant, OHS Branch, Clayton campus –Original to be filed in area.**